



MANUFACTURING COMPANY LIMITED

BLACKHORSE ROAD, LETCHWORTH, HERTS, SG6 1HB, ENGLAND
Tel: + 44 (0) 1462 686221 www.stoddard.co.uk Fax: + 44 (0) 1462 480711

PRODUCT SAFETY DATA SHEET

Product: COMPO Silicone Polishers, RA Shanks
Classification: Class 2a Medical Device

Schedule of Products (Ref):

S151, S152, S153, S154	POLISHERS, COMPO
------------------------	------------------

Description: Stoddard abrasive COMPO Silicon polishers are an Interlaced Silicon, end polymerised product joined to metal shank.

Component

Abrasive Compo polishers (working part) head.

Description: Interlaced Silicon, end polymerised with abrasive grains and coloured pigments.

Composition:

Silicon mix 35%
Abrasive grains 60%
Colour pigments, 5%

The materials contained in this product are non-hazardous or non-hazardous at the levels of inclusion. This material may also contain trace elements, which are non-hazardous or non-hazardous at the levels of inclusion.

Physical data:

Physical condition: Solid Substance
Colour: Various
Odour: N/A
Change of Condition: N/A
Density: Approx. 1.3 – 1.5 g/cm³ (20 °C)
Vapour Pressure: N/A
Viscosity: Solid Substance
Solubility in water: Insoluble
pH value: N/A
Flash point: N/A
Ignition temperature: Not specified
Explosion Limits: N/A

Dangerous Materials category: No dangerous Substance.

First Aid Measures:

Contact with eyes: In case of contact with eyes, rinse eyes immediately with large amounts of water (10 to 15 Minutes). Call a Doctor if necessary.
Contact with Skin: Wash thoroughly with Water.
Ingestion: N/A, insoluble in Stomach and Intestines.
Inhalation: N/A

Product Safety Data Sheets

First fighting procedures:

Extinguishing Agents: Water, Carbon Dioxide, Foam, Powder

Special measures in case of Unintended release:

None Required.

Handling and Storage:

Handling: Use of this product is only recommended by a suitably qualified dental surgeon or Hygienist
Hygiene: Take customary hygiene measures of the trade.
Storage: Store in a clean and dry place. Keep away from sources of heat. Do not expose to direct sunlight. Always store in genuine packaging.

Exploitation limitation and Protective equipment:

Technical measures: Good general ventilation is deemed sufficient.
Limit value monitoring: No limit value specified.

Personal Protective equipment:

Respirator: Respirator (Dust)
Hand Protection: None
Eye Protection: Protective goggles or face shield.

Stability and Reactivity:

Thermal Disintegration: None at designated use.
Hazardous disintegration Products: None under normal storage or use
Dangerous reactions: None
Further Information: Protect against exposure to light during storage.

Toxicology Information:

Acute Toxicity: None
Chronic Toxicity: None

Ecology Information:

Water Imperilment category: 0 (No water imperilment)

Disposal:

Dispose as domestic waste

Transport:

Overland transport: GGVE / GGVS: No Classification.
RID / ADR: N/A
(Follow normal rules on overland transportation of hazardous substances)
Sea transport: GGVSee: No Classification.
(Follow normal rules on Ocean transportation of hazardous substances)
Air Transport: ICAO / IATA-DGR: No Classification.
(Follow normal rules on Air transportation of hazardous substances)

Further Information:

The products are designed for use in direct contact with Human beings and are subject to EC rule 93/42 for medical devices.

Product Safety Data Sheets

DISCLAIMER OF LIABILITY

The information in this PSDS was obtained from sources, which we believe are reliable. However, the information is provided without any warranty, express or implied, regarding its correctness.

The conditions or methods of handling, storage, use or disposal of the product are beyond our control and may be beyond our knowledge.

For this and other reasons, we do not assume responsibility and expressly disclaim liability for loss, damage or expense arising out of or in any way connected with the handling, storage, use or disposal of the product.

This PSDS was prepared and is to be used only for this product. If the product is used as a component in another product, this PSDS information may not be applicable.

Emergency Phone No: (0044) 1462 686221